

## **Manifold Molecular Iodine Cells**

(Product ID: I2M-5, I2M-10)

To create a sealed starved cell with a flexible set point, a glass cell with an attached cold-finger and vacuum port is constructed. The vacuum port and cold-finger include stopcocks. The cell is evacuated and cold-finger filled with Iodine is brought to the desired vapor pressure (cold-finger operating temperature). The stem between the cold-finger and cell is then closed by closing the stopcock, isolating the Iodine in the cell body and fixing the number density. The cell is then operated 10-20 °C above the cold-finger set temperature and the Iodine in the cell is a super-heated vapor with a set number density. The result is a molecular cell with a very stable absorption spectra.

Iodine-vapor cells with a manifold are 3-in.-dia, 5-in.-long or 10-in.-long Pyrex cells. Standard cells are manufactured with 1/4 gm of iodine (actual pressure determined by temperature of water cooling jacket at the time of operation). In addition, transitions can be pressure broadened with buffer gas supplied by user through a fill-port. Both models are EAR99.




Packaged I2M-10 with thermocouples and heater.

## **SPECIFICATIONS FOR 12M-5**

Length	125-mm
Diameter	76-mm
Housing	Anodized Aluminum
Mounting Options	
Maximum Operating Temperature	130 °C
Thermocouple	Type T
Set Point	

## SPECIFICATIONS FOR 12M-10

Length	250-mm
Diameter	
Housing	Anodized Aluminum
Mounting Options	1/4 -20 Threaded
Maximum Operating Temperature	130 °C
Thermocouple	Type T
Set Point	



Transmission versus frequency for 5-inch iodine cell at a range  $of\ cold\hbox{-}finger\ (vapor\ pressure)\ set\ points.$ 



Side view of I2M-10 cell design.

Innovative Scientific Solutions, Inc. 7610 McEwen Road, Dayton, OH 45459 Phone (937) 630-3012 Fax (937) 630-3015

email: spalluconi@innssi.com

web: www.psp-tsp.com